


ARTHRO**VIMS**<sup>®</sup>

# MOSAIC PLASTY



## DESCRIPTION

Mosaic plasty is an osteochondral autograft. Osteochondral grafts of various diameters are harvested from unstressed cartilage knee areas and then re-implanted to reconstruct the articular surface of a lesional area. The ArthroVims mosaic plasty system is ideally suited to indications involving small areas of cartilage damage in knee and ankle.



### Simple, ergonomic system

Different diameters available (6mm/8mm/10mm) to cover all surgical indications

## PRODUCTS & INSTRUMENTATIONS

System suitable for arthroscopic or mini-open surgery.

Designation	Reference
<b>Single-use sampling/insertion kit</b>	
System plot donor/recipient Ø 6mm	MZPK-0006
System plot donor/recipient Ø 8mm	MZPK-0008
System plot donor/recipient Ø 10mm	MZPK-0010
<b>Restorable ancillary</b>	
Handle donor	HNDD- 001
Handle recipient	HNDR-0001
Tamp	TAMP-0001
Recipient Reamer Ø 6mm	RPRM-0006
Recipient Reamer Ø 8mm	RPRM-0008
Recipient Reamer Ø 10mm	RPRM-0010
Insert (2 pieces in set)	DRIN-0001
Wrench for handle	WRNC-0001
Hammer	HMMR-0001
Container	MSPS-0001