



ALL SUTURE ANCHOR

SOFT ANCHOR

SA

ARTHRO **VIMS**[®]

SOFT ANCHOR





SHOULDER
ROTATOR CUFF,
BICEPS TENODESIS AND
REINSERTION

Ø2,9
mm



SHOULDER
SLAP, BANKART
INSTABILITY



ANKLE
TENDON AND LIGAMENT
REINSERTION

Ø1,8
mm



ANKLE
TENDON AND LIGAMENT
REINSERTION



SHOULDER
SLAP, BANKART
INSTABILITY

Ø1,5
mm



HAND
TENDON AND LIGAMENT
REINSERTION

Ø1,2
mm



SUTURE ADVANTAGES

ARTHROVIMS SOFT ANCHOR are associated with 2 or 3 USP2 sutures
(made of very high molecular weight polyethylene UHMWPE)

- HIGH TENSILE STRENGTH
- EASY REALIZATION OF KNOTS
- STRENGTH AND SAFETY OF SUTURES



STORAGE

ARTHROVIMS® SOFT ANCHOR must be stored in its original unopened packaging, in a clean and dry place, at ambient temperature.

PRODUCTS & INSTRUMENTATION

ARTHROVIMS® SOFT ANCHOR is packaged in a double pouch and sterilized by ethylene oxide. Single use.

Do not re-sterilize.

Instrumentation for arthroscopic procedures is available.

ARTHROVIMS® SOFT ANCHOR must be used with the corresponding instrumentation supplied by VIMS. For any further information, please refer to the IFU.

DESIGNATION	SIZE	REFERENCE
ARTHROVIMS SOFT ANCHOR	2,9	VIMS-029
ARTHROVIMS SOFT ANCHOR with needles	2,9 needle	VIMS-029N
ARTHROVIMS SOFT ANCHOR	1,8	VIMS-018
ARTHROVIMS SOFT ANCHOR with needles	1,8 needle	VIMS-018N
ARTHROVIMS SOFT ANCHOR	1,5	VIMS-015
ARTHROVIMS SOFT ANCHOR with needles	1,5 needle	VIMS-015N
ARTHROVIMS SOFT ANCHOR with needles	1,2 needle	VIMS-012N

INSTRUMENTATION

ARTHROVIMS® SOFT ANCHOR must be used with the specific instrumentation provided by VIMS.



CE 1984
CLASS IIb / ANCHOR

CE
CLASS I / INSTRUMENTATION